

A 19  
Lid  
For fixing delicate wash items.



- Perforated plate 7 x 7 x 3 mm
- Suitable for insert E 969

EAN: 4002513298758 / Material number: 05746210 / Old Material Number: 69796901D

<b>Product category</b>	
Large-capacity washer-disinfectors, laboratory	•
<b>Machine affiliation</b>	
PG 8527	•
PLW 8615	•
PLW 8616	•
PLW 8617	•
<b>Product category</b>	
Lid	•
<b>Product properties</b>	
Material	Stainless steel
Colour	Stainless Steel
<b>Area of application</b>	
Reprocessing of laboratory glassware	•
<b>Standard electrical connection</b>	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
<b>Dimensions and weight</b>	
External dimensions, net height in mm	18
External dimensions, net width in mm	251
External dimensions, net depth in mm	351
External dimensions, gross height in mm	70
External dimensions, gross width in mm	275
External dimensions, gross depth in mm	515
Net weight in kg	0.55
Gross weight in kg	0.70