

E 757

Mobile unit

for the optimum loading of large volume laboratory glassware.

EAN: 4002513298772 / Material number: 05746290 / Old Material Number: 69575701D



- Equipped with 6 nozzles
- Equipped with a height-adjustable holding frame
- Maximum vertical clearance from star overlay: 610 mm mm

Product category	
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
G 7825	•
Product category	
Other mobile units	•
Module for large-volume laboratory glassware	•
Product properties	
Material	Stainless steel, Plastic
Colour	Stainless Steel, Black
Maximum vertical clearance in mm	610
Area of application	
Reprocessing of laboratory glassware	•
Capacity	
Large-volume laboratory glassware [number]	6
Standard electrical connection	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
Dimensions and weight	
External dimensions, net height in mm	346
External dimensions, net width in mm	530
External dimensions, net depth in mm	600
External dimensions, gross height in mm	675
External dimensions, gross width in mm	565
External dimensions, gross depth in mm	650
Net weight in kg	10.90
Gross weight in kg	10.90