

E 957

Mobile unit

for the optimum loading of up to 12 large volume laboratory glassware.



- Equipped with 12 nozzles with Ø 10 mm
- Height adjustable holding frame
- Magnetic strip for automatic mobile unit recognition
- Max. vertical clearance from star overlay: 615 mm mm

EAN: 4002513298789 / Material number: 05746300 / Old Material Number: 69595701D

<b>Product category</b>	
Large-capacity washer-disinfectors, laboratory	•
<b>Machine affiliation</b>	
PG 8527	•
<b>Product category</b>	
Other mobile units	•
Module for large-volume laboratory glassware	•
<b>Product properties</b>	
Material	Stainless steel, Plastic
Colour	Stainless Steel, Black
Maximum vertical clearance in mm	615
<b>Area of application</b>	
Reprocessing of laboratory glassware	•
<b>Capacity</b>	
Large-volume laboratory glassware [number]	12
<b>Standard electrical connection</b>	
Number of phases	0
Voltage in V	0-0
Frequency in Hz	0-0
Total rated load in kW	0.00-0.00
Fuse rating in A	0-0
<b>Dimensions and weight</b>	
External dimensions, net height in mm	353
External dimensions, net width in mm	640
External dimensions, net depth in mm	790
External dimensions, gross height in mm	390
External dimensions, gross width in mm	665
External dimensions, gross depth in mm	805
Net weight in kg	16.98
Gross weight in kg	17.45