

E 131/1

Insert

for optimum loading of 5 mesh and kidney trays.



- Suitable for 5 mesh or kidney trays
- Equipped with 6 holders, height: 160 mm mm, distance: 80 mm mm

EAN: 4002513317336 / Material number: 05771030 / Old Material Number: 69513102D

Product category	
Washer-disinfectors, medical	•
Machine affiliation	
PG 8535	•
PG 8536	•
PG 8581	•
PG 8562	•
PG 8582	•
PG 8582 CD	•
PG 8591	•
PG 8592	•
Product category	
Insert for mesh trays	•
Insert for kidney dishes/basins	•
Product properties	
Material	Stainless steel
Colour	Stainless Steel
Distance between rails in mm	80
Rail height in mm	160
Area of application	
Reprocessing of kidney dishes	•
Reprocessing of surgical instruments	•
Capacity	
Kidney dishes [number]	5
Standard electrical connection	
Number of phases	0
Voltage in V	0-0
Frequency in Hz	0-0
Total rated load in kW	0.00-0.00
Fuse rating in A	0-0
Dimensions and weight	
External dimensions, net height in mm	168
External dimensions, net width in mm	180
External dimensions, net depth in mm	495
External dimensions, gross height in mm	185
External dimensions, gross width in mm	200
External dimensions, gross depth in mm	530
Net weight in kg	0.61
Gross weight in kg	0.90