


E 965  
 Insert  
 for the optimum loading of beakers 250 - 600 ml.



- Equipped with 15 x 3 stay bars for beakers from 250 ml to 600 ml

EAN: 4002513369229 / Material number: 05848290 / Old Material Number: 69596501D

<b>Product category</b>	
Large-capacity washer-disinfectors, laboratory	•
<b>Machine affiliation</b>	
PG 8527	•
<b>Product category</b>	
Insert/module for beakers	•
<b>Product properties</b>	
Material	Stainless steel
Colour	Stainless Steel
<b>Area of application</b>	
Reprocessing of laboratory glassware	•
<b>Capacity</b>	
Beakers [quantity]	98
<b>Standard electrical connection</b>	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
<b>Dimensions and weight</b>	
External dimensions, net height in mm	173
External dimensions, net width in mm	357
External dimensions, net depth in mm	522
External dimensions, gross height in mm	280
External dimensions, gross width in mm	540
External dimensions, gross depth in mm	560
Net weight in kg	2.07
Gross weight in kg	2.86