


E 960/1

Insert

for the optimum loading of wide necked Erlenmeyer flasks.



- Equipped with 20 large and 26 small spring hooks
- Suitable for Erlenmeyer flasks & measuring cylinders

EAN: 4002513400953 / Material number: 05892360 / Old Material Number: 69596002D

Product category	
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
PG 8527	•
Product category	
Insert/module for wide-necked glassware	•
Product properties	
Material	Stainless steel
Colour	Stainless Steel
Area of application	
Reprocessing of laboratory glassware	•
Capacity	
Erlenmeyer flasks [number]	46
Standard electrical connection	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
Dimensions and weight	
External dimensions, net height in mm	185
External dimensions, net width in mm	357
External dimensions, net depth in mm	522
External dimensions, gross height in mm	280
External dimensions, gross width in mm	530
External dimensions, gross depth in mm	560
Net weight in kg	0.80
Gross weight in kg	0.90