


O 177/1

Upper basket

for optimum cleaning with 26 nozzles and holding of various inserts.



- Right side for different inserts
- Left side with 26 nozzles and 7 funnels supplied loose
- Equipped with height adjustable frame
- Vertical clearance 230/205 mm
- Height adjustable: – 20 / – 40 mm

EAN: 4002513412413 / Material number: 05913050 / Old Material Number: 69117702D

<b>Product category</b>	
Thermal disinfectors	•
Washer-disinfectors, medical	•
<b>Machine affiliation</b>	
PG 8535	•
PG 8536	•
<b>Product category</b>	
Carrier/basket carrier	•
Upper basket	•
<b>Product properties</b>	
Material	Stainless steel, Plastic
Colour	Stainless Steel, Black
<b>Area of application</b>	
Reprocessing of ENT instruments	•
Reprocessing of dental instruments	•
<b>Capacity</b>	
Transmission instruments [number]	11
Lumened instruments [number]	22
<b>Standard electrical connection</b>	
Number of phases	0
Voltage in V	0-0
Frequency in Hz	0-0
Total rated load in kW	0.00-0.00
Fuse rating in A	0-0
<b>Dimensions and weight</b>	
External dimensions, net height in mm	263
External dimensions, net width in mm	498
External dimensions, net depth in mm	455
External dimensions, gross height in mm	280
External dimensions, gross width in mm	545
External dimensions, gross depth in mm	560
Net weight in kg	3.70
Gross weight in kg	4.52