

E 441/1

Insert

for the optimum loading of Micro instruments.



- The mesh width on the base is 1.7 mm
- The sides are closed, stackable
- Sections w/ 6 adjustable dividers for positioning of instruments
- Can be steam sterilised at 121 °C/134 °C

EAN: 4002513460476 / Material number: 05988240 / Old Material Number: 69544102D

<b>Product category</b>	
Thermal disinfectors	•
Washer-disinfectors, medical	•
<b>Machine affiliation</b>	
G 7831	•
PG 8535	•
PG 8536	•
PG 8581	•
PG 8562	•
PG 8582	•
PG 8582 CD	•
PG 8591	•
PG 8592	•
<b>Product category</b>	
Insert/module for MIS instruments	•
<b>Product properties</b>	
Material	Stainless steel, Plastic
Colour	Stainless Steel, Black
<b>Area of application</b>	
Reprocessing of MIS instruments	•
<b>Standard electrical connection</b>	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
<b>Dimensions and weight</b>	
External dimensions, net height in mm	60
External dimensions, net width in mm	183
External dimensions, net depth in mm	284
External dimensions, gross height in mm	100
External dimensions, gross width in mm	210
External dimensions, gross depth in mm	315
Net weight in kg	0.77
Gross weight in kg	0.92