

FP 20

Filter plates

for the reprocessing of contra angles and hand pieces when using AUF 1/AUF 2.

EAN: 4002513505498 / Material number: 06059480 / Old Material Number: 69751201D



- To avoid blockages by micro particles
- Porosity 2
- 20 filter plates per pack

Product category	
Thermal disinfectors	•
Large-capacity washer-disinfectors, medical	•
Machine affiliation	
G 7823	•
G 7824	•
G 7825	•
G 7826	•
PG 8527	•
PG 8528	•
G 7831	•
PG 8535	•
PG 8536	•
Product properties	
Colour	White
Standard electrical connection	
Number of phases	0
Voltage in V	0-0
Frequency in Hz	0-0
Total rated load in kW	0.00-0.00
Fuse rating in A	0-0
Dimensions and weight	
External dimensions, net height in mm	30
External dimensions, net width in mm	100
External dimensions, net depth in mm	125
External dimensions, gross height in mm	28
External dimensions, gross width in mm	125
External dimensions, gross depth in mm	167
Diameter in mm	30
Net weight in kg	0.06
Gross weight in kg	0.10