


E 745/1
Module
for the optimum loading of up to 104 pipettes.



- Suitable for 104 pipettes up to 540 mm
- Holding frame compartment size 16 x 16 mm

EAN: 4002513735406 / Material number: 06233580 / Old Material Number: 69574502D

Product category	
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
G 7825	•
Product category	
Insert/module for pipettes	•
Product properties	
Material	Stainless steel, Plastic
Colour	Stainless Steel, Black
Area of application	
Reprocessing of pipettes/phials	•
Capacity	
Pipettes [number]	104
Standard electrical connection	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
Dimensions and weight	
External dimensions, net height in mm	288
External dimensions, net width in mm	492
External dimensions, net depth in mm	496
External dimensions, gross height in mm	415
External dimensions, gross width in mm	530
External dimensions, gross depth in mm	590
Net weight in kg	24.60
Gross weight in kg	25.00