

E 529

Mobile unit

for optimum loading of 4-5 ophthalmic theatre sets.



- 2 levels for mesh trays
- 1 level with 16 or 30 injector nozzles
- Connection option for E 478/1 to the inlet pipe
- Equipped with connection for hot air drying
- Magnetic strip for mobile unit recognition (without magnets ML)

EAN: 4002513792553 / Material number: 06415950 / Old Material Number: 69552901D

Product category	
Large-capacity washer-disinfectors, medical	•
Machine affiliation	
G 7823	•
G 7824	•
Product category	
Mobile unit for ophthalmology instruments	•
Product properties	
Material	Stainless steel
Colour	Stainless Steel
Area of application	
Reprocessing of ophthalmology instruments	•
Capacity	
Theatre sets for micro instruments [number]	4-5
Lumened instruments [number]	46
Standard electrical connection	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
Dimensions and weight	
External dimensions, net height in mm	471
External dimensions, net width in mm	530
External dimensions, net depth in mm	645
External dimensions, gross height in mm	550
External dimensions, gross width in mm	640
External dimensions, gross depth in mm	650
Net weight in kg	14.00
Gross weight in kg	21.00