


U 184/1

Lower basket

for the optimum loading of centrifuge tubes, phials and test tubes.



- Equipped with 96 injector nozzles
- Only suitable for use with upper baskets O 175 / O 187 / O 184
- Vertical clearance: Approx. 170 mm

EAN: 4002513805659 / Material number: 06439180 / Old Material Number: 69118403D

<b>Product category</b>	
Washer-disinfectors, laboratory	•
<b>Machine affiliation</b>	
PG 8536	•
<b>Product category</b>	
Injector basket	•
Lower basket	•
<b>Product properties</b>	
Material	Stainless steel, Plastic
Colour	Stainless Steel, Black
Maximum vertical clearance in mm	170
<b>Area of application</b>	
Reprocessing of pipettes/phials	•
Reprocessing of test tubes	•
<b>Capacity</b>	
Test tubes [number]	96
Phials [number]	96
<b>Standard electrical connection</b>	
Number of phases	0
Voltage in V	0-0
Frequency in Hz	0-0
Total rated load in kW	0.00-0.00
Fuse rating in A	0-0
<b>Dimensions and weight</b>	
External dimensions, net height in mm	250
External dimensions, net width in mm	515
External dimensions, net depth in mm	531
External dimensions, gross height in mm	270
External dimensions, gross width in mm	525
External dimensions, gross depth in mm	540
Net weight in kg	5.00
Gross weight in kg	6.00
<b>Companion pack</b>	
ID 90, injector nozzle, Ø 2.5 mm, length 90 mm [number]	96