

E 436/3

Mobile unit

for optimum loading of approx. 3 anaesthesia sets.



- Suitable for about 3 anaesthesia sets.
- 8 nozzles with spring holders for breathing tubes up to 1.5 m
- Storage on spiral
- Equipped with connection for hot air drying
- Bracket for ML/2 for mobile unit recognition / spray arm control

EAN: 4002513805673 / Material number: 06439240 / Old Material Number: 69543604D

<b>Product category</b>	
Washer-disinfectors, medical	•
<b>Machine affiliation</b>	
PG 8536	•
<b>Product category</b>	
Mobile unit for anaesthetic instruments	•
<b>Product properties</b>	
Material	Stainless steel
Colour	Stainless Steel
<b>Area of application</b>	
Reprocessing of anaesthetic instruments	•
<b>Capacity</b>	
AN sets [number]	3
<b>Standard electrical connection</b>	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
<b>Dimensions and weight</b>	
External dimensions, net height in mm	507
External dimensions, net width in mm	535
External dimensions, net depth in mm	535
External dimensions, gross height in mm	510
External dimensions, gross width in mm	540
External dimensions, gross depth in mm	520
Net weight in kg	6.60
Gross weight in kg	7.00
<b>Companion pack</b>	
E 430/1 1/3 mesh tray [number]	1
E 431 injector nozzle for concertina hose, Ø 8 x 193 mm [number]	1
E 432 holder for concertina breathing tubes [number]	1
E 433 holder for silicone breathing tubes [number]	2
E 434 holder for paediatric breathing tubes [number]	1
E 466 injector nozzle 8 x 333 mm, for breathing bags [number]	6
E 496 injector nozzle 4 x 120 mm [number]	10