

E 427

Insert

for the optimum loading of 6 laryngoscopes.



- Suitable for 6 laryngoscopes

EAN: 4002513808209 / Material number: 06442120 / Old Material Number: 69542701D

<b>Product category</b>	
Washer-disinfectors, medical	•
Large-capacity washer-disinfectors, medical	•
<b>Machine affiliation</b>	
G 7823	•
G 7824	•
G 7825	•
G 7826	•
PG 8528	•
PG 8535	•
PG 8536	•
PG 8562	•
PG 8582	•
PG 8582 CD	•
PG 8592	•
<b>Product category</b>	
Insert/module for anaesthetic instruments	•
<b>Product properties</b>	
Material	Stainless steel
Colour	Stainless Steel
<b>Area of application</b>	
Reprocessing of anaesthetic instruments	•
<b>Standard electrical connection</b>	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
<b>Dimensions and weight</b>	
External dimensions, net height in mm	92
External dimensions, net width in mm	210
External dimensions, net depth in mm	134
External dimensions, gross height in mm	164
External dimensions, gross width in mm	245
External dimensions, gross depth in mm	250
Net weight in kg	0.90
Gross weight in kg	1.00