

E 479

Holder

for micro instruments up to 4 mm Ø.



- For use in trays with 5 mm mesh
- 50 items per bag

EAN: 4002513918175 / Material number: 06553750 / Old Material Number: 69747901D

| | |
|---|-------------|
| Product category | |
| Washer-disinfectors, medical | • |
| Large-capacity washer-disinfectors, medical | • |
| Machine affiliation | |
| G 7823 | • |
| G 7824 | • |
| G 7825 | • |
| G 7826 | • |
| PG 8527 | • |
| PG 8528 | • |
| PG 8535 | • |
| PG 8536 | • |
| PG 8562 | • |
| PG 8582 | • |
| PG 8582 CD | • |
| PG 8592 | • |
| Product category | |
| Holder | • |
| Product properties | |
| Material | Plastic |
| Colour | Transparent |
| Area of application | |
| Reprocessing of MIS instruments | • |
| Reprocessing of ophthalmology instruments | • |
| Standard electrical connection | |
| Number of phases | 0 |
| Voltage in V | 0-0 |
| Frequency in Hz | 0-0 |
| Total rated load in kW | 0.00-0.00 |
| Fuse rating in A | 0-0 |
| Dimensions and weight | |
| External dimensions, net height in mm | 5 |
| External dimensions, net width in mm | 21 |
| External dimensions, net depth in mm | 21 |
| External dimensions, gross height in mm | 95 |
| External dimensions, gross width in mm | 110 |
| External dimensions, gross depth in mm | 80 |
| Net weight in kg | 0.19 |
| Gross weight in kg | 0.23 |