


A 11/1

Insert

For the gentle cleaning of delicate items by interrupting the spray jet.



- With perforations 7 x 7 mm
- For use in upper and lower basket
- Recommended for use in PG 8536

EAN: 4002513975680 / Material number: 06629660 / Old Material Number: 69501102D

Product category	
Washer-disinfectors, laboratory	•
Machine affiliation	
PG 8535	•
PG 8536	•
PG 8583	•
PG 8583 CD	•
PG 8593	•
Product category	
Retaining frame	•
Product properties	
Material	Stainless steel
Colour	Stainless Steel
Area of application	
Reprocessing of laboratory glassware	•
Standard electrical connection	
Number of phases	0
Voltage in V	0-0
Frequency in Hz	0-0
Total rated load in kW	0.00-0.00
Fuse rating in A	0-0
Dimensions and weight	
External dimensions, net height in mm	1
External dimensions, net width in mm	433
External dimensions, net depth in mm	433
External dimensions, gross height in mm	16
External dimensions, gross width in mm	470
External dimensions, gross depth in mm	470
Net weight in kg	0.80
Gross weight in kg	1.00