

E 430/1

Insert

For optimum accommodation of different instruments.



- Comprising wire mesh with 5 mm mesh width

EAN: 4002513984422 / Material number: 06643200 / Old Material Number: 69743003D

Product category	
Thermal disinfectors	•
Washer-disinfectors, medical	•
Machine affiliation	
G 7831	•
PG 8535	•
PG 8536	•
PG 8581	•
PG 8562	•
PG 8582	•
PG 8582 CD	•
PG 8591	•
PG 8592	•
Product category	
Mesh tray	•
Product properties	
Material	Stainless steel
Colour	Stainless Steel
Area of application	
Reprocessing of surgical instruments	•
Standard electrical connection	
Number of phases	0
Voltage in V	0-0
Frequency in Hz	0-0
Total rated load in kW	0.00-0.00
Fuse rating in A	0-0
Dimensions and weight	
External dimensions, net height in mm	40
External dimensions, net width in mm	150
External dimensions, net depth in mm	445
External dimensions, gross height in mm	100
External dimensions, gross width in mm	160
External dimensions, gross depth in mm	540
Net weight in kg	0.36
Gross weight in kg	0.48