

E 405/1

Mobile unit

for optimum loading of up to 38 pipettes, connection for hot air drying.



- Suitable for 38 pipettes in 3 rows
- 1st Row: 10 pipettes x 100 ml (up to 550 mm), bar distance 20 mm
- 2nd Row: 14 pipettes x 25 ml, bar distance 26 mm
- 3rd Row: 14 pipettes x 10 ml, bar distance 26 mm
- Bracket for ML/2 for mobile unit recognition / spray arm control

EAN: 4002514087078 / Material number: 06727920 / Old Material Number: 69540502D

Product category	
Washer-disinfectors, laboratory	•
Machine affiliation	
PG 8535	•
PG 8536	•
Product category	
Mobile unit with connection for hot-air drying	•
Mobile unit for pipettes	•
Product properties	
Material	Stainless steel, Plastic
Colour	Stainless Steel, Black
Number of rows	3
Area of application	
Reprocessing of pipettes/phials	•
Capacity	
Pipettes [number]	38
Standard electrical connection	
Number of phases	0
Voltage in V	0-0
Frequency in Hz	0-0
Total rated load in kW	0.00-0.00
Fuse rating in A	0-0
Dimensions and weight	
External dimensions, net height in mm	502
External dimensions, net width in mm	535
External dimensions, net depth in mm	517
External dimensions, gross height in mm	500
External dimensions, gross width in mm	540
External dimensions, gross depth in mm	530
Net weight in kg	6.20
Gross weight in kg	8.10