

U 874/2

Lower basket

for standard DIN mesh trays and Miele inserts.



- Vertical clearance can vary with different upper baskets

EAN: 4002514110967 / Material number: 06756110 / Old Material Number: 69187404D

Product category	
Washer-disinfectors, laboratory	•
Thermal disinfectors	•
Washer-disinfectors, medical	•
Machine affiliation	
PG 8535	•
PG 8536	•
Product category	
Carrier/basket carrier	•
Lower basket	•
Product properties	
Material	Stainless steel
Colour	Stainless Steel
Area of application	
Reprocessing of laboratory glassware	•
Reprocessing of MIS instruments	•
Reprocessing of GYN instruments	•
Reprocessing of ENT instruments	•
Reprocessing of dental instruments	•
Reprocessing of theatre shoes	•
Reprocessing of basins	•
Reprocessing of kidney dishes	•
Standard electrical connection	
Number of phases	0
Dimensions and weight	
External dimensions, net height in mm	50
External dimensions, net width in mm	534
External dimensions, net depth in mm	515
External dimensions, gross height in mm	125
External dimensions, gross width in mm	555
External dimensions, gross depth in mm	555
Net weight in kg	2.50
Gross weight in kg	3.00