


U 874/2

Lower basket

for standard DIN mesh trays and Miele inserts.



- Vertical clearance can vary with different upper baskets

EAN: 4002514110967 / Material number: 06756110 / Old Material Number: 69187404D

<b>Product category</b>	
Washer-disinfectors, laboratory	•
Thermal disinfectors	•
Washer-disinfectors, medical	•
<b>Machine affiliation</b>	
PG 8535	•
PG 8536	•
<b>Product category</b>	
Carrier/basket carrier	•
Lower basket	•
<b>Product properties</b>	
Material	Stainless steel
Colour	Stainless Steel
<b>Area of application</b>	
Reprocessing of laboratory glassware	•
Reprocessing of MIS instruments	•
Reprocessing of GYN instruments	•
Reprocessing of ENT instruments	•
Reprocessing of dental instruments	•
Reprocessing of theatre shoes	•
Reprocessing of basins	•
Reprocessing of kidney dishes	•
<b>Standard electrical connection</b>	
Number of phases	0
<b>Dimensions and weight</b>	
External dimensions, net height in mm	50
External dimensions, net width in mm	534
External dimensions, net depth in mm	515
External dimensions, gross height in mm	125
External dimensions, gross width in mm	555
External dimensions, gross depth in mm	555
Net weight in kg	2.50
Gross weight in kg	3.00