

KWT 4154 UG-1



EAN: 4002514149226 / Material number: 06804170 /
Old Material Number: 36415403GB

General information	
Model identifier	KWT 4154 UG-1
Manufacturer material number	06804170
EAN number	4002514149226
Marketing class	Platinum
Appliance category	
Wine conditioning unit	•
Construction type	
Undercounter appliance	•
Door hinge	right
Convertible door hinging	•
Design	
Casing colour	Stainless Steel
Front colour	Glass door/handle
Wine compartment lighting	LED illumination
User convenience	
DynaCool	•
Active AirClean Filter	•
Ice/Water dispenser with mains water connection	No
Number of FlexiFrame racks	3
Control system	
Operating concept	ComfortControl
Fridge compartment can be switched off	No
No. of temperature zones	2
Min. temperature setting range in °C	5
Safety	
Acoustic door alarm	•
Acoustic temperature alarm	•
Technical data	
Niche width in mm min.	600
Niche width in mm max.	600
Niche height in mm min.	820
Niche height in mm max.	870
Niche depth in mm	580
Appliance width in mm	597
Appliance height in mm	817
Appliance depth in mm	575
Weight in kg	52.0
Climate range	SN-ST
Wine compartment in l	126
Of which pantry compartment in l	126
0.75 l Bordeaux bottle capacity (no. of bottles)	38 Flaschen
Storage time in case of fault in h	0
Current consumption (mA)	800
Voltage in V	220-240
Fuse rating in A	13
Number of phases	1
Frequency in Hz	50

KWT 4154 UG-1




EAN: 4002514149226 / Material number: 06804170 /
Old Material Number: 36415403GB

Length of supply lead in m

2.0

KWT 4154 UG-1



KWT4154UG, KWT4154UG-1 (Sketch)

Footnote: , *If the neighbouring kitchen unit doors are taller than the wine cooler door, the wine cooler can be aligned to the unit doors using a filler panel.