

E 139/1

Insert

For the optimum loading of test tubes up to 12 x 200 mm.



- Suitable for test tubes up to 12 x 200 mm
- Divided into 6 surfaces
- Mesh width 9 x 9 mm
- For use in lower basket
- Including lid A 13

EAN: 4002514231082 / Material number: 06907660 / Old Material Number: 69513902D

Product category	
Washer-disinfectors, laboratory	•
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
G 7825	•
PG 8527	•
PG 8535	•
PG 8536	•
PG 8504	•
PG 8583	•
PG 8583 CD	•
PG 8593	•
PLW 8615	•
PLW 8616	•
PLW 8617	•
PLW 6011	•
PLW 6111	•
PLW 7111	•
PLW 8505	•
Product category	
Insert/module for test tubes	•
Product properties	
Material	Stainless steel
Colour	Stainless Steel
Number of compartments	6
Area of application	
Reprocessing of test tubes	•
Reprocessing of laboratory glassware	•
Capacity	
Test tubes [number]	200
Standard electrical connection	
Number of phases	0
Voltage in V	0-0
Frequency in Hz	0-0
Total rated load in kW	0.00-0.00
Fuse rating in A	0-0
Dimensions and weight	
External dimensions, net height in mm	243
External dimensions, net width in mm	200
External dimensions, net depth in mm	230
External dimensions, gross height in mm	245
External dimensions, gross width in mm	245
External dimensions, gross depth in mm	265
Net weight in kg	1.26
Gross weight in kg	1.29
Companion pack	
A 13 lid [number]	1