

PG 8597

Twin-tank water softener

Twin-tank water softener for the continuous removal of soft water up to 40 °dH

EAN: 4002514256283 / Material number: 06943120 / Old Material Number: 69859701D



- Free-standing appliance with rollers, plastic outer casing
- Output: continuous 19 l/min, max. 30 l/min
- 2 containers with 4.5 l resin each, 1 container for 20 kg salt
- Freestanding appliance with rollers, top loaded
- continuous 19 l/min, max. 30 l/min
- 1 container for 20 kg salt, 2 containers for 4.5 l resin
- Volume controlled two-chamber system, current-less operation
- Volume controlled two-chamber system, current-less operation
- The appliance is equipped with a water connection

Product category	
Tank dishwashers	•
Washer-disinfectors, laboratory	•
Large-capacity washer-disinfectors, laboratory	•
Thermal disinfectors	•
Washer-disinfectors, medical	•
Large-capacity washer-disinfectors, medical	•
Machine affiliation	
G 8066	•
PG 8067	•
PG 8164	•
PG 8165	•
PG 8166	•
PG 8169	•
PG 8172	•
PG 8172 Eco	•
PTD 901	•
PTD 701	•
PTD 702	•
PTD 703	•
G 7823	•
G 7824	•
G 7825	•
G 7826	•
PG 8527	•
PTD 704	•
PG 8528	•
G 7831	•
PG 8535	•
PG 8536	•
Product category	
Water treatment	•
Product properties	
Material	Plastic
Colour	White
Standard electrical connection	
Number of phases	0
Voltage in V	0-0
Frequency in Hz	0-0
Total rated load in kW	0.00-0.00
Fuse rating in A	0-0
Dimensions and weight	
External dimensions, net height in mm	570
External dimensions, net width in mm	360
External dimensions, net depth in mm	360
External dimensions, gross height in mm	680
External dimensions, gross width in mm	375
External dimensions, gross depth in mm	380
Net weight in kg	19.00
Gross weight in kg	26.70