

E 905/1

Module

for optimum loading of MIS instruments, such as arthroscopy and urology.



- Holds up to 16 hollow instruments
- Subdivision for arthroscopy and laparoscopy, for example.

EAN: 4002514256306 / Material number: 06943240 / Old Material Number: 69590502D

Product category	
Washer-disinfectors, medical	•
Large-capacity washer-disinfectors, medical	•
Machine affiliation	
G 7823	•
G 7824	•
G 7825	•
G 7826	•
PG 8527	•
PG 8528	•
PG 8536	•
Product category	
Insert/module for MIS instruments	•
Product properties	
Material	Stainless steel
Colour	Stainless Steel
Area of application	
Reprocessing of MIS instruments	•
Capacity	
Lumened instruments [number]	16
Standard electrical connection	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
Dimensions and weight	
External dimensions, net height in mm	40
External dimensions, net width in mm	461
External dimensions, net depth in mm	510
External dimensions, gross height in mm	170
External dimensions, gross width in mm	540
External dimensions, gross depth in mm	560
Net weight in kg	3.01
Gross weight in kg	3.62
Companion pack	
E 336 MIBO irrigation sleeve, Ø 11 mm, length 121 mm [number]	1
E 362 blind stopper [number]	2
E 442 irrigation sleeve, Ø 11 mm, length 121 mm, with cap Ø 6 mm [number]	1
E 447 female adapter, for male Luer Lock [number]	2
E 448 silicone hose with female Luer Lock adapter [number]	4
E 449 male adapter, for female Luer Lock [number]	2
E 452 injector nozzle 2.5 x 60 mm [number]	4
E 453 injector nozzle 4 x 110 mm with retaining clip [number]	3
E 454 injector nozzle for trocar sleeves 10–15 mm [number]	3
E 464 holder for injector nozzle E 454 [number]	1
E 472 spring clip for injector nozzles Ø 4 mm [number]	1

E 905/1

Module

for optimum loading of MIS instruments, such as arthroscopy and urology.

EAN: 4002514256306 / Material number: 06943240 / Old Material Number: 69590502D



E 907/1 mesh tray with lid for small parts [number]

1