

E 912/2

Mobile unit

For optimum loading of 4 ISO containers including lids.



- Container size: 180 x 400 x 600 mm
- Magnetic strip for mobile unit recognition (incl. magnets ML)

EAN: 4002514257976 / Material number: 06943860 / Old Material Number: 69591203D

<b>Product category</b>	
Large-capacity washer-disinfectors, medical	•
<b>Machine affiliation</b>	
PG 8527	•
PG 8528	•
<b>Product category</b>	
Mobile unit for containers	•
<b>Product properties</b>	
Material	Stainless steel
Colour	Stainless Steel
<b>Area of application</b>	
Reprocessing of containers	•
<b>Capacity</b>	
ISO containers height 18 cm [number]	4
<b>Standard electrical connection</b>	
Number of phases	0
Voltage in V	0-0
Frequency in Hz	0-0
Total rated load in kW	0.00-0.00
Fuse rating in A	0-0
<b>Dimensions and weight</b>	
External dimensions, net height in mm	455
External dimensions, net width in mm	640
External dimensions, net depth in mm	790
External dimensions, gross height in mm	650
External dimensions, gross width in mm	650
External dimensions, gross depth in mm	800
Net weight in kg	13.00
Gross weight in kg	22.00
<b>Companion pack</b>	
Magnetic strip for automatic mobile unit recognition (incl. magnets) [number]	1