

E 915/3

Mobile unit

for optimum loading of approx. 7 AN-sets or 20 breathing tubes and accessories.



- Holds breathing tubes on a spiral rack
- Injector system
- Equipped with connection for hot air drying
- Magnetic strip for mobile unit recognition (incl. magnets)

EAN: 4002514257983 / Material number: 06943870 / Old Material Number: 69591504D

<b>Product category</b>	
Large-capacity washer-disinfectors, medical	•
<b>Machine affiliation</b>	
PG 8527	•
PG 8528	•
<b>Product category</b>	
Mobile unit with connection for hot-air drying	•
Mobile unit for anaesthetic instruments	•
<b>Product properties</b>	
Material	Stainless steel
Colour	Stainless Steel
<b>Area of application</b>	
Reprocessing of anaesthetic instruments	•
<b>Capacity</b>	
AN sets [number]	7
<b>Standard electrical connection</b>	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
<b>Dimensions and weight</b>	
External dimensions, net height in mm	637
External dimensions, net width in mm	640
External dimensions, net depth in mm	790
External dimensions, gross height in mm	780
External dimensions, gross width in mm	670
External dimensions, gross depth in mm	810
Net weight in kg	25.00
Gross weight in kg	25.00
<b>Companion pack</b>	
E 142 insert 1/2, DIN mesh tray [number]	1
E 431 injector nozzle for concertina hose, Ø 8 x 193 mm [number]	1
E 466 injector nozzle 8 x 333 mm, for breathing bags [number]	10
E 496 injector nozzle 4 x 120 mm [number]	18
E 916 holder for 5 concertina breathing tubes [number]	1
E 917 holder for 5 silicone breathing tubes [number]	4
E 918 holder for 5 paediatric breathing tubes [number]	1
Spring holders, Ø 12 mm, for breathing tubes up to 1.5 m long [number]	20
Magnetic strip for automatic mobile unit recognition (incl. magnets) [number]	1