


O 188/2  
Upper basket  
for the optimum loading of different inserts.



- Vertical clearance: 165 mm, height adjustable +/- 20 mm
- With integrated magnetic spray arm
- Magnet for use with spray arm sensors for PG 8535, PG 8536

EAN: 4002514556178 / Material number: 07333030 / Old Material Number: 69118802D

<b>Product category</b>	
Washer-disinfectors, laboratory	•
Washer-disinfectors, medical	•
<b>Machine affiliation</b>	
PG 8535	•
PG 8536	•
<b>Product category</b>	
Carrier/basket carrier	•
Upper basket	•
<b>Product properties</b>	
Material	Stainless steel
Colour	Stainless Steel
<b>Standard electrical connection</b>	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
<b>Dimensions and weight</b>	
External dimensions, net height in mm	280
External dimensions, net width in mm	215
External dimensions, net depth in mm	475
External dimensions, gross height in mm	285
External dimensions, gross width in mm	540
External dimensions, gross depth in mm	555
Net weight in kg	2.40
Gross weight in kg	3.16