

O 190/2

Upper basket

for the optimum loading of different inserts.



- Vertical clearance: 215 mm mm
- With integrated magnetic spray arm
- Magnet for use with spray arm sensors for PG 8535, PG 8536

EAN: 4002514556185 / Material number: 07333050 / Old Material Number: 69119002D

Product category	
Washer-disinfectors, laboratory	•
Washer-disinfectors, medical	•
Machine affiliation	
PG 8535	•
PG 8536	•
Product category	
Carrier/basket carrier	•
Upper basket	•
Product properties	
Material	Stainless steel
Colour	Stainless Steel
Maximum vertical clearance in mm	215
Standard electrical connection	
Number of phases	0
Voltage in V	0-0
Frequency in Hz	0-0
Total rated load in kW	0.00-0.00
Fuse rating in A	0-0
Dimensions and weight	
External dimensions, net height in mm	265
External dimensions, net width in mm	531
External dimensions, net depth in mm	475
External dimensions, gross height in mm	285
External dimensions, gross width in mm	540
External dimensions, gross depth in mm	555
Net weight in kg	2.40
Gross weight in kg	3.19