


O 167/1

Upper basket

for the optimum loading of insoles for theatre shoes.



- Suitable for up to 28 insoles
- With integrated magnetic spray arm
- Magnet for use with spray arm sensors for PG 8535, PG 8536

EAN: 4002514556215 / Material number: 07333080 / Old Material Number: 69116703D

<b>Product category</b>	
Washer-disinfectors, medical	•
<b>Machine affiliation</b>	
PG 8535	•
PG 8536	•
<b>Product category</b>	
Insert for theatre shoes	•
Upper basket	•
<b>Product properties</b>	
Material	Stainless steel
Colour	Stainless Steel
<b>Area of application</b>	
Reprocessing of theatre shoes	•
<b>Standard electrical connection</b>	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
<b>Dimensions and weight</b>	
External dimensions, net height in mm	195
External dimensions, net width in mm	480
External dimensions, net depth in mm	475
External dimensions, gross height in mm	210
External dimensions, gross width in mm	495
External dimensions, gross depth in mm	540
Net weight in kg	2.40
Gross weight in kg	3.23