

O 173/1

Upper basket

for the optimum loading of theatre shoes.



- Suitable for 8 theatre shoes up to size 41
- With integrated magnetic spray arm
- Magnet for use with spray arm sensors for PG 8535, PG 8536

EAN: 4002514556222 / Material number: 07333090 / Old Material Number: 69117302D

Product category	
Washer-disinfectors, medical	•
Machine affiliation	
PG 8535	•
PG 8536	•
Product category	
Insert for theatre shoes	•
Upper basket	•
Product properties	
Material	Stainless steel
Colour	Stainless Steel
Area of application	
Reprocessing of theatre shoes	•
Capacity	
Theatre shoes [number]	8
Standard electrical connection	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
Dimensions and weight	
External dimensions, net height in mm	195
External dimensions, net width in mm	531
External dimensions, net depth in mm	475
External dimensions, gross height in mm	540
External dimensions, gross width in mm	270
External dimensions, gross depth in mm	525
Net weight in kg	2.40
Gross weight in kg	3.30