

MF 27/28-1

Transfer trolley

for loading & unloading washer-&disinfectors.



- Equipped with 4 castors, 2 of which can be locked
- Simple and convenient handling thanks to rounded handles
- Dockable or lockable to rolling conveyor etc. with both sides
- Infeed height 850 mm, -100 mm, +150 mm
- With removable drip tray H 70, W 603, D 866 mm

EAN: 4002514594026 / Material number: 07397640 / Old Material Number: 69200111D

| | |
|--|-----------------|
| Product category | |
| Large-capacity washer-disinfectors, laboratory | • |
| Large-capacity washer-disinfectors, medical | • |
| Machine affiliation | |
| PG 8527 | • |
| PG 8528 | • |
| Product category | |
| Storage/transport system | • |
| Product properties | |
| Material | Stainless steel |
| Colour | Stainless Steel |
| Standard electrical connection | |
| Number of phases | 0 |
| Voltage in V | 0-0 |
| Frequency in Hz | 0-0 |
| Total rated load in kW | 0.00-0.00 |
| Fuse rating in A | 0-0 |
| Dimensions and weight | |
| External dimensions, net height in mm | 1050 |
| External dimensions, net width in mm | 740 |
| External dimensions, net depth in mm | 930 |
| External dimensions, gross height in mm | 1010 |
| External dimensions, gross width in mm | 730 |
| External dimensions, gross depth in mm | 950 |
| Net weight in kg | 32.00 |
| Gross weight in kg | 45.00 |