

E 943/2

Module

for the optimum loading of narrow necked glassware 100 - 500 ml.



- 32 nozzles (E 351) 4 x 160 mm with clips (E 353)

EAN: 4002514605043 / Material number: 07459400 / Old Material Number: 69594303D

Product category	
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
PG 8527	•
PG 8528	•
Product category	
Insert/module for narrow-necked glassware	•
Product properties	
Material	Stainless steel
Colour	Stainless Steel
Area of application	
Reprocessing of laboratory glassware	•
Capacity	
Narrow-necked glassware [number]	32
Standard electrical connection	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
Dimensions and weight	
External dimensions, net height in mm	190
External dimensions, net width in mm	558
External dimensions, net depth in mm	352
External dimensions, gross height in mm	180
External dimensions, gross width in mm	560
External dimensions, gross depth in mm	550
Net weight in kg	6.50
Gross weight in kg	7.50