

E 947/2

Module

for cleaning 88 centrifuge tubes, phials or test tubes, for example.



- Equipped with 88 nozzles, 2.5 x 110 mm

EAN: 4002514605074 / Material number: 07459430 / Old Material Number: 69594703D

Product category	
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
PG 8527	•
PG 8528	•
Product category	
Insert/module for test tubes	•
Product properties	
Material	Stainless steel, Plastic
Colour	Stainless Steel, Black
Area of application	
Reprocessing of pipettes/phials	•
Capacity	
Test tubes [number]	88
Standard electrical connection	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
Dimensions and weight	
External dimensions, net height in mm	170
External dimensions, net width in mm	558
External dimensions, net depth in mm	352
External dimensions, gross height in mm	290
External dimensions, gross width in mm	530
External dimensions, gross depth in mm	560
Net weight in kg	5.00
Gross weight in kg	6.00