

DOS S 20  
Bellows dispenser pump  
For dispensing liquid acid media.



- For dispensing rinsing agent, neutralising agent
- For retro-fitting by Miele Service

EAN: 4002514641904 / Material number: 07514600 /  
Old Material Number: 69226004D

<b>Product category</b>	
Laboratory washers	•
Washer-disinfectors, medical	•
<b>Machine affiliation</b>	
PG 8536	•
<b>Product category</b>	
Dispenser system	•
<b>Standard electrical connection</b>	
Number of phases	0
Voltage in V	0-0
Frequency in Hz	0-0
Total rated load in kW	0.00-0.00
Fuse rating in A	0-0
<b>Dimensions and weight</b>	
External dimensions, gross height in mm	245
External dimensions, gross width in mm	258
External dimensions, gross depth in mm	270
Net weight in kg	2.45
Gross weight in kg	2.91