

E 935/2

Mobile unit

for the optimum cleaning of instruments on 3 levels.



- Loading capacity: 9 DIN mesh trays*
- Equipped with 2 integrated spray arms
- The loading areas for the individual levels are different
- Magnetic strip for mobile unit recognition (incl. magnets)

* DIN mesh tray 485 x 254 x 50 mm

EAN: 4002514816432 / Material number: 07765780 / Old Material Number: 69593503D

Product category	
Large-capacity washer-disinfectors, laboratory	•
Large-capacity washer-disinfectors, medical	•
Machine affiliation	
PG 8527	•
PG 8528	•
Product category	
Mobile unit for inserts	•
Mobile unit for mesh trays	•
Product properties	
Material	Stainless steel
Colour	Stainless Steel
Area of application	
Reprocessing of laboratory glassware	•
Reprocessing of surgical instruments	•
Capacity	
DIN mesh trays [number]	9
SPRI mesh trays [number]	6
Standard electrical connection	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
Dimensions and weight	
External dimensions, net height in mm	524
External dimensions, net width in mm	640
External dimensions, net depth in mm	790
External dimensions, gross height in mm	720
External dimensions, gross width in mm	830
External dimensions, gross depth in mm	780
Net weight in kg	26.50
Gross weight in kg	30.00
Companion pack	
Magnetic strip for automatic mobile unit recognition (incl. magnets) [number]	1