

E 975/2

Mobile unit

for the optimum cleaning of instruments on 2 levels.



- Equipped with an integrated spray arm
- Magnetic strip for automatic mobile unit recognition
- The vertical clearance for both levels is different

EAN: 4002514816449 / Material number: 07765790 / Old Material Number: 69597503D

<b>Product category</b>	
Large-capacity washer-disinfectors, laboratory	•
<b>Machine affiliation</b>	
PG 8527	•
PG 8528	•
<b>Product category</b>	
Mobile unit for inserts	•
Mobile unit for mesh trays	•
<b>Product properties</b>	
Material	Stainless steel
Colour	Stainless Steel
<b>Area of application</b>	
Reprocessing of laboratory glassware	•
Reprocessing of surgical instruments	•
<b>Capacity</b>	
DIN mesh trays [number]	6
SPRI mesh trays [number]	4
<b>Standard electrical connection</b>	
Number of phases	0
Voltage in V	0-0
Frequency in Hz	0-0
Total rated load in kW	0.00-0.00
Fuse rating in A	0-0
<b>Dimensions and weight</b>	
External dimensions, net height in mm	426
External dimensions, net width in mm	636
External dimensions, net depth in mm	790
External dimensions, gross height in mm	650
External dimensions, gross width in mm	550
External dimensions, gross depth in mm	640
Net weight in kg	10.80
Gross weight in kg	12.80
<b>Companion pack</b>	
Magnetic strip for automatic mobile unit recognition (incl. magnets) [number]	1