

APH 407
Network cable
for connection to the network, length 5 m.



- Cable category 6
- Colour: Transparent silver

EAN: 4002514940588 / Material number: 07951480 / Old Material Number: 68815705D

Product category	
Washer-disinfectors, laboratory	•
Large-capacity washer-disinfectors, laboratory	•
Thermal disinfectors	•
Washer-disinfectors, medical	•
Large-capacity washer-disinfectors, medical	•
Small sterilisers	•
Large sterilisers	•
Machine affiliation	
G 7823	•
G 7824	•
G 7825	•
G 7826	•
PG 8527	•
PG 8528	•
G 7831	•
PG 8535	•
PG 8536	•
PG 8581	•
PG 8582	•
PG 8582 CD	•
PG 8583	•
PG 8583 CD	•
PG 8591	•
PG 8592	•
PG 8593	•
PS 1201 B	•
PS 1201 B EXCELLENCE	•
PS 1202 B	•
PS 5662 V	•
PS 5666 V	•
PS 5669 V	•
PS 5962 H	•
PS 5966 H	•
PS 5968 H	•
PS 5969 H	•
Product category	
Process documentation/operating documentation/ communication	•
Product properties	
Colour	Silver
Area of application	
Reprocessing of pipettes/phials	•
Reprocessing of laboratory glassware	•
Reprocessing of slides	•
Reprocessing of butyrometers	•
Reprocessing of pharmacy/infusion bottles	•
Reprocessing of anaesthetic instruments	•

APH 407
 Network cable
 for connection to the network, length 5 m.



EAN: 4002514940588 / Material number: 07951480 / Old Material Number: 68815705D

Reprocessing of MIS instruments	•
Reprocessing of ophthalmology instruments	•
Reprocessing of robotic instruments	•
Reprocessing of GYN instruments	•
Reprocessing of ENT instruments	•
Reprocessing of dental instruments	•
Reprocessing of theatre shoes	•
Reprocessing of basins	•
Reprocessing of kidney dishes	•
Reprocessing of baby bottles	•
Reprocessing of surgical instruments	•
Reprocessing of lumened instruments	•
Reprocessing of orthopaedic surgery instruments	•
Reprocessing of intramedullary nails/bone marrow drills	•
Reprocessing of trays	•
Reprocessing of mesh trays	•
Reprocessing of measuring cylinders	•
Reprocessing of containers, e.g. soda lime containers	•
Reprocessing of different utensils	•
Dimensions and weight	
External dimensions, gross height in mm	12
External dimensions, gross width in mm	100
External dimensions, gross depth in mm	150
Net weight in kg	0.13
Gross weight in kg	0.13