

E 428

Mobile unit

for the optimum loading of instruments used in robot-assisted surgery.



- Suitable for 6 instruments for robot-assisted surgery
- Has dual connection and 1 reserve connection per instrument
- The ROBOTVARIO programme is required for the use of the mobile unit\*

\* The chemicals are not included and must be ordered separately.

EAN: 4002514948805 / Material number: 07973140 / Old Material Number: 69542801D

<b>Product category</b>	
Washer-disinfectors, medical	•
<b>Machine affiliation</b>	
PG 8536	•
<b>Product category</b>	
Mobile unit for robotic instruments	•
<b>Product properties</b>	
Material	Stainless steel, Plastic
Colour	Stainless Steel, Black, Transparent
<b>Area of application</b>	
Reprocessing of robotic instruments	•
<b>Capacity</b>	
Robotics instruments [number]	6
<b>Dimensions and weight</b>	
External dimensions, net height in mm	518
External dimensions, net width in mm	533
External dimensions, net depth in mm	513
External dimensions, gross height in mm	530
External dimensions, gross width in mm	550
External dimensions, gross depth in mm	550
Net weight in kg	2.40
Gross weight in kg	3.30