


E 197

Insert

For optimum accommodation of different instruments.



- Mesh width on the base 3 mm, sides 1.7 mm, lid 8 mm
- Can be steam sterilised at 121 °C/134 °C

EAN: 4002514953380 / Material number: 07984090 / Old Material Number: 68519701D

| | |
|---|-----------------|
| Product classification | |
| Thermal disinfectors | • |
| Machine affiliation | |
| G 7831 | • |
| PG 8581 | • |
| PG 8591 | • |
| PWD 8531 | • |
| PWD 8532 | • |
| Product category | |
| Mesh tray | • |
| Product properties | |
| Material | Stainless steel |
| Colour | Stainless Steel |
| Area of application | |
| Reprocessing of dental instruments | • |
| Standard electrical connection | |
| Number of phases | 0 |
| Voltage in V | 0 |
| Frequency in Hz | 0 |
| Total rated load in kW | 0.00 |
| Fuse rating in A | 0 |
| Dimensions and weight | |
| External dimensions, net height in mm | 42 |
| External dimensions, net width in mm | 150 |
| External dimensions, net depth in mm | 225 |
| External dimensions, gross height in mm | 75 |
| External dimensions, gross width in mm | 170 |
| External dimensions, gross depth in mm | 240 |
| Net weight in kg | 0.55 |
| Gross weight in kg | 0.62 |