


E 807

Insert

for optimum loading of 3 mesh or kidney trays.



- 4 holders (3 compartments), H 165 mm, spacing 68 mm mm

EAN: 4002514951638 / Material number: 07984220 / Old Material Number: 69580701D

<b>Product classification</b>	
Thermal disinfectors	•
Washer-disinfectors, medical	•
<b>Machine affiliation</b>	
G 7831	•
PG 8581	•
PG 8591	•
PWD 8531	•
PWD 8532	•
PWD 8534	•
<b>Product category</b>	
Insert for mesh trays	•
Insert for kidney dishes/basins	•
<b>Product properties</b>	
Material	Stainless steel
Colour	Stainless Steel
Number of compartments	3
Distance between rails in mm	68
<b>Area of application</b>	
Reprocessing of dental instruments	•
<b>Capacity</b>	
Diverse mesh trays [number]	3
Kidney dishes [number]	3
<b>Standard electrical connection</b>	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
<b>Dimensions and weight</b>	
External dimensions, net height in mm	144
External dimensions, net width in mm	140
External dimensions, net depth in mm	290
External dimensions, gross height in mm	480
External dimensions, gross width in mm	320
External dimensions, gross depth in mm	130
Net weight in kg	0.31
Gross weight in kg	0.47