

APH 330

Adapter gender changer

to connect two same-pole serial connections.



- Compatible for use with HAWO heat-seal units
- 9-pin plug, m/m

EAN: 4002514954219 / Material number: 07986870 / Old Material Number: 68815506D

Product category	
Process documentation/operating documentation/communication	•
Product properties	
Colour	Silver
Area of application	
Reprocessing of anaesthetic instruments	•
Reprocessing of MIS instruments	•
Reprocessing of ophthalmology instruments	•
Reprocessing of robotic instruments	•
Reprocessing of ENT instruments	•
Reprocessing of dental instruments	•
Reprocessing of kidney dishes	•
Reprocessing of surgical instruments	•
Reprocessing of lumened instruments	•
Reprocessing of orthopaedic surgery instruments	•
Reprocessing of intramedullary nails/bone marrow drills	•
Dimensions and weight	
External dimensions, net height in mm	0
External dimensions, net width in mm	0
External dimensions, net depth in mm	0
Net weight in kg	0.02
Gross weight in kg	0.02