

E 339/1

Insert

for the optimum loading of up to 12 tray dividers and trays.



- Suitable for 12 tray dividers or trays
- Equipped with 13 holders (12 compartments).
- Maximum tray size 290 x 20 mm

EAN: 4002515094860 / Material number: 09233840 / Old Material Number: 69533902D

<b>Product category</b>	
Thermal disinfectors	•
Washer-disinfectors, medical	•
<b>Machine affiliation</b>	
G 7831	•
PG 8535	•
PG 8536	•
PG 8581	•
PG 8562	•
PG 8582	•
PG 8582 CD	•
PG 8591	•
PG 8592	•
<b>Product category</b>	
Insert for trays	•
<b>Product properties</b>	
Material	Stainless steel
Colour	Stainless Steel
Number of compartments	12
<b>Capacity</b>	
Trays [quantity]	12
<b>Standard electrical connection</b>	
Number of phases	0
Voltage in V	0-0
Frequency in Hz	0-0
Total rated load in kW	0.00-0.00
Fuse rating in A	0-0
<b>Dimensions and weight</b>	
External dimensions, net height in mm	115
External dimensions, net width in mm	305
External dimensions, net depth in mm	498
External dimensions, gross height in mm	480
External dimensions, gross width in mm	320
External dimensions, gross depth in mm	130
Net weight in kg	0.80
Gross weight in kg	1.22