

APH 331

Connection cable

For connection of Melag Vacuklav 24, 31 and Autoclav 25.



- In addition to connection cable type 1
- Length 3 m
- For connecting the machine to Segosoft Miele Edition

EAN: 4002515350973 / Material number: 09573240 / Old Material Number: 68812505D

Product category	
Small sterilisers	•
Product category	
Process documentation/operating documentation/communication	•
Product properties	
Colour	Grey
Area of application	
Reprocessing of pharmacy/infusion bottles	•
Reprocessing of anaesthetic instruments	•
Reprocessing of MIS instruments	•
Reprocessing of ophthalmology instruments	•
Reprocessing of GYN instruments	•
Reprocessing of dental instruments	•
Reprocessing of surgical instruments	•
Reprocessing of lumened instruments	•
Reprocessing of orthopaedic surgery instruments	•
Reprocessing of mesh trays	•
Reprocessing of containers, e.g. soda lime containers	•
Reprocessing of different utensils	•
Dimensions and weight	
External dimensions, net height in mm	0
External dimensions, net width in mm	0
External dimensions, net depth in mm	0
External dimensions, gross height in mm	12
External dimensions, gross width in mm	300
External dimensions, gross depth in mm	400
Net weight in kg	0.12
Gross weight in kg	0.13