

B 2826 FashionMaster Steam ironing system



EAN: 4002515434666 / Material number: 09819810 /
Old Material Number: 13282621GB

General information	
Model identifier	B 2826 FashionMaster
Manufacturer material number	09819810
EAN number	4002515434666
Design	
Iron colour	Anthracite/lotus white
Ironing board cover colour	Anthracite
Perfect ironing results	
Honeycomb soleplate	•
Steam pressure in bar	4.0
Steam quantity in g/min	100
Vertical steam	•
Steam starting from temperature level	2
Active ironing table (inflation and suction function)	•
Suction function	•
2-stage fan	•
Non-stick soleplate	•
Time-saving ironing	
Water container volume in l	1.25
Removable water container	•
Permanently refillable water container	•
Multi-layer ironing	•
User convenience	
All-in-one	•
1-2 Lift system	•
Variable height adjustment with gas pressure stay, in mm	830 - 1020
Castors gentle on the floor	•
Ironing board with Comfort Zone	•
Fan safety switch-on/switch-off	•
Fan memory function	•
Automatic descaling	•
Automatic rinsing	•
Storage compartment for iron	•
Storage compartment for mains cable and accessories	•
Safety	
CoolDown function	•
AutoOff function	•
Technical data	
Dimensions in mm (width)	470
Dimensions in mm (height)	970
Dimensions in mm (depth)	1480
Dimensions in mm (width), folded	470
Dimensions in mm (height), folded	1280
Dimensions in mm (depth), folded	370
Ironing surface area in mm	1200 x 400
Frequency in Hz	50 - 60
Accessories included	
Non-stick soleplate	•

B 2826 FashionMaster Steam ironing system



EAN: 4002515434666 / Material number: 09819810 /
Old Material Number: 13282621GB

Premium ironing table cover

Steam hose holder

•
•