


A 101

Upper basket

Height adjustable, for the optimum loading of inserts.



- Vertical clearance 160 +/- 30 mm, width 475 mm, depth 443 mm
- With integrated spray arm

EAN: 4002515460627 / Material number: 09862340 / Old Material Number: 69610100D

Product category	
Fresh-water dishwashers	•
Washer-disinfectors, laboratory	•
Washer-disinfectors, medical	•
Machine affiliation	
PG 8063	•
PG 8562	•
PG 8582	•
PG 8582 CD	•
PG 8504	•
PG 8583	•
PG 8583 CD	•
PG 8592	•
PG 8593	•
Product category	
Upper basket	•
Carrier/basket carrier	•
Upper basket	•
Product properties	
Material	Stainless steel
Colour	Stainless Steel
Area of application	
Reprocessing of laboratory glassware	•
Reprocessing of GYN instruments	•
Reprocessing of ENT instruments	•
Reprocessing of basins	•
Reprocessing of kidney dishes	•
Reprocessing of surgical instruments	•
Reprocessing of different utensils	•
Reprocessing of breathing air regulators	•
Standard electrical connection	
Number of phases	0
Voltage in V	0-0
Frequency in Hz	0-0
Total rated load in kW	0.00-0.00
Fuse rating in A	0-0
Dimensions and weight	
External dimensions, net height in mm	206
External dimensions, net width in mm	528
External dimensions, net depth in mm	527
External dimensions, gross height in mm	230
External dimensions, gross width in mm	555
External dimensions, gross depth in mm	580
Net weight in kg	3.00
Gross weight in kg	3.79