


A 104

Upper basket

For the optimum loading of instruments for MIS, arthroscopy and urology.



- Right side for holding 11 hollow instruments
- Additional Luer lock connectors for hollow instruments
- With spray arm: Mesh trays can be placed on the lower carrier

EAN: 4002515460658 / Material number: 09862380 / Old Material Number: 69610400D

| | |
|---|-----------------|
| Product category | |
| Washer-disinfectors, medical | • |
| Machine affiliation | |
| PG 8562 | • |
| PG 8582 | • |
| PG 8582 CD | • |
| PG 8592 | • |
| Product category | |
| Basket for MIS instruments | • |
| Injector basket | • |
| Upper basket | • |
| Product properties | |
| Material | Stainless steel |
| Colour | Stainless Steel |
| Area of application | |
| Reprocessing of MIS instruments | • |
| Capacity | |
| DIN mesh trays [number] | 1 |
| Lumened instruments [number] | 11 |
| Standard electrical connection | |
| Number of phases | 0 |
| Voltage in V | 0-0 |
| Frequency in Hz | 0-0 |
| Total rated load in kW | 0.00-0.00 |
| Fuse rating in A | 0-0 |
| Dimensions and weight | |
| External dimensions, net height in mm | 332 |
| External dimensions, net width in mm | 528 |
| External dimensions, net depth in mm | 526 |
| External dimensions, gross height in mm | 270 |
| External dimensions, gross width in mm | 525 |
| External dimensions, gross depth in mm | 540 |
| Net weight in kg | 4.00 |
| Gross weight in kg | 6.50 |