


A 104

Upper basket

For the optimum loading of instruments for MIS, arthroscopy and urology.



- Right side for holding 11 hollow instruments
- Additional Luer lock connectors for hollow instruments
- With spray arm: Mesh trays can be placed on the lower carrier

EAN: 4002515460658 / Material number: 09862380 / Old Material Number: 69610400D

<b>Product category</b>	
Washer-disinfectors, medical	•
<b>Machine affiliation</b>	
PG 8562	•
PG 8582	•
PG 8582 CD	•
PG 8592	•
<b>Product category</b>	
Basket for MIS instruments	•
Injector basket	•
Upper basket	•
<b>Product properties</b>	
Material	Stainless steel
Colour	Stainless Steel
<b>Area of application</b>	
Reprocessing of MIS instruments	•
<b>Capacity</b>	
DIN mesh trays [number]	1
Lumened instruments [number]	11
<b>Standard electrical connection</b>	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
<b>Dimensions and weight</b>	
External dimensions, net height in mm	332
External dimensions, net width in mm	528
External dimensions, net depth in mm	526
External dimensions, gross height in mm	270
External dimensions, gross width in mm	525
External dimensions, gross depth in mm	540
Net weight in kg	4.00
Gross weight in kg	6.50