

A 150

Lower basket

For the optimum loading of injector modules or inserts.



- Lower basket with two module connections
- for holding up to 2 injector modules or inserts
- Automatic closing of the connections when not in use
- Clearances: width 494, length 440 mm

EAN: 4002515460672 / Material number: 09862400 / Old Material Number: 69615000D

<b>Product category</b>	
Washer-disinfectors, laboratory	•
<b>Machine affiliation</b>	
PG 8504	•
PG 8583	•
PG 8583 CD	•
PG 8593	•
<b>Product category</b>	
Basket with module connection	•
Lower basket	•
<b>Product properties</b>	
Material	Stainless steel, Plastic
Colour	Stainless Steel, Black
<b>Area of application</b>	
Reprocessing of pipettes/phials	•
Reprocessing of laboratory glassware	•
Reprocessing of butyrometers	•
Reprocessing of measuring cylinders	•
<b>Standard electrical connection</b>	
Number of phases	0
Voltage in V	0-0
Frequency in Hz	0-0
Total rated load in kW	0.00-0.00
Fuse rating in A	0-0
<b>Dimensions and weight</b>	
External dimensions, net height in mm	154
External dimensions, net width in mm	529
External dimensions, net depth in mm	546
External dimensions, gross height in mm	175
External dimensions, gross width in mm	545
External dimensions, gross depth in mm	565
Net weight in kg	2.61
Gross weight in kg	3.29