

A 202  
Mobile unit  
for the optimum loading of 4 DIN mesh trays.



- Suitable for 4 DIN mesh trays on 2 levels
- Insert frames: Removable for the reprocessing of tall wash items
- 6 additional connection options via Luer lock connection
- With integrated spray arm
- Loading volumes for the different levels differ from each other

EAN: 4002515460719 / Material number: 09862440 / Old Material Number: 69620200D

<b>Product category</b>	
Washer-disinfectors, laboratory	•
Washer-disinfectors, medical	•
<b>Machine affiliation</b>	
PG 8562	•
PG 8582	•
PG 8582 CD	•
PG 8583	•
PG 8583 CD	•
PG 8592	•
PG 8593	•
<b>Product category</b>	
Mobile unit for inserts	•
Mobile unit for mesh trays	•
<b>Product properties</b>	
Material	Stainless steel, Plastic
Colour	Stainless Steel, Black
<b>Area of application</b>	
Reprocessing of laboratory glassware	•
Reprocessing of surgical instruments	•
<b>Capacity</b>	
DIN mesh trays [number]	4
<b>Standard electrical connection</b>	
Number of phases	0
Voltage in V	0-0
Frequency in Hz	0-0
Total rated load in kW	0.00-0.00
Fuse rating in A	0-0
<b>Dimensions and weight</b>	
External dimensions, net height in mm	223
External dimensions, net width in mm	529
External dimensions, net depth in mm	542
External dimensions, gross height in mm	235
External dimensions, gross width in mm	550
External dimensions, gross depth in mm	560
Net weight in kg	5.50
Gross weight in kg	6.74