

A 203

Mobile unit

for the optimum loading of MIS instruments.



- Suitable for holding and for E 451, E 457, E 460, E 473, E 444

EAN: 4002515460726 / Material number: 09862450 / Old Material Number: 69620300D

<b>Product category</b>	
Washer-disinfectors, medical	•
<b>Machine affiliation</b>	
PG 8562	•
PG 8582	•
PG 8582 CD	•
PG 8592	•
<b>Product category</b>	
Mobile unit for MIS instruments	•
<b>Product properties</b>	
Material	Stainless steel
Colour	Stainless Steel
<b>Area of application</b>	
Reprocessing of MIS instruments	•
<b>Standard electrical connection</b>	
Number of phases	0
Voltage in V	0-0
Frequency in Hz	0-0
Total rated load in kW	0.00-0.00
Fuse rating in A	0-0
<b>Dimensions and weight</b>	
External dimensions, net height in mm	520
External dimensions, net width in mm	529
External dimensions, net depth in mm	542
External dimensions, gross height in mm	270
External dimensions, gross width in mm	525
External dimensions, gross depth in mm	540
Net weight in kg	9.40
Gross weight in kg	11.00
<b>Companion pack</b>	
E 336 MIBO irrigation sleeve, Ø 11 mm, length 121 mm [number]	3
E 362 blind stopper [number]	2
E 442 irrigation sleeve, Ø 11 mm, length 121 mm, with cap Ø 6 mm [number]	15
E 443 irrigation sleeve, Ø 11 mm, length 121 mm, with cap Ø 10 mm [number]	5
E 445 – 12 covers for E 442 irrigation sleeve, opening with Ø 6 mm [number]	1
E 446 – 12 covers for E 443 irrigation sleeve, opening with Ø 10 mm [number]	1
E 447 female adapter, for male Luer Lock [number]	3
E 448 silicone hose with female Luer Lock adapter [number]	6
E 449 male adapter, for female Luer Lock [number]	5
E 451 insert 1/6 mesh tray [number]	1
E 452 injector nozzle 2.5 x 60 mm [number]	3
E 453 injector nozzle 4 x 110 mm with retaining clip [number]	8
E 454 injector nozzle for trocar sleeves 10–15 mm [number]	6
E 456 spring stay for MIS instruments [number]	4
E 464 holder for injector nozzle E 454 [number]	3